

OXYGEN CYLINDER DURATION TIME

Oxygen Cylinder Duration Chart

D Cylinder	
LBS Pressure	10 LPM
2,000	39 min
1,500	28 min
1,000	18 min
500	7 min
E Cylinder	
LBS Pressure	10 LPM
2,000	63 min
1,500	45 min
1,000	28 min
500	11 min
M Cylinder	
LBS Pressure	10LPM
2,000	5 hr 30 min
1,500	4 hr
1,000	2 hr 30 min
500	1 hr

*Please note: All times shown are approximate.
Actual time is dependent on factors such as temperature, elevation, etc.*



D Size cylinder
Capacity 416 Liters of O₂



E Size cylinder
Capacity 682 Liters of O₂



M Size cylinder
Capacity 3450 Liters of O₂